



MöllerMedical Aesthetics Solutions

Surgest
Medical



Workstation Pro plus

Five reasons

1 TLA infiltration, liposuction and fat re-infiltration - **Workstation Pro plus** optimally supports you.

2 The included **Liposat® Pro plus** in combination with **Vibrasat® Pro** handles fat in the gentlest way - Because viable fat matters!

3 The unique **Anti-Hypothermia Solution** – **Thermia Pro** increases the patient safety during the surgery.

4 Enjoy a super silent OR and be able to fully concentrate on your patient - Maximum reduction of the noise level.

5 The integrated 3-Pedal Footswitch allows you to perform your procedure exactly as you need it - Offer your patients the best possible results.

Consists of



Vibrasat® Pro
REF 00003920



Liposat® Pro plus
REF 00003974
incl. 1 Pedal-Footswitch



Vacusat® power
REF 00002250



Thermia Pro mounted
REF 00004253



3-Pedal Footswitch | 5m wire
REF 00003981



Overflow protection
REF 00002299



Hydrophobic filter
REF 00002297

REF: 00004272

Are you shaping the right way? - The **Workstation Pro plus**, offers everything you need to achieve outstanding results and revolutionize your way of shaping.

With the combination of **Vibrasat® Pro** and **Liposat® Pro plus**, caring about the valuable fat cells becomes as easy as the fat re-infiltration routine itself.

Key benefits at one glance

Best with prior TLA infiltration

Rollable Easy operation

Designed for
**Liposuction and
fat re-infiltration**

Included foot pedals for hand-free operation

Overflow protection **Perfectly matched devices**

For small and large volume

Noise optimisation Space-saving and handy design

Hydrophobic filter Slims the workflow

Easy cleanability Reduces preparing time



Surgest
Medical

Distribuidor Exclusivo en España
Teléfono de atención al cliente: 935 895 350
surgest.com



MöllerMedical Aesthetics Solutions

Liposat® Pro plus

Surgest
Medical



Five reasons

1 Gently and accurately handle fat cells. To maximize the viability of fat cells the influence of the pump and the tube set on the viability of the fat cells is reduced to a minimum.

2 According to an in-house in-vitro study, **Liposat® Pro plus** achieves a 92 % or more survival rate for infiltrated fat cells.

3 The versatile two-in-one pump is developed for both TLA and fat infiltration.

4 Ensure reliable results and a high level of safety for your patients, maximizing their satisfaction with the procedure.

5 **Liposat® Pro plus** takes care of the volume management, so you can focus exclusively on your patient.



REF: 00003974

Liposat® Pro plus - because we care about viable fat cells!



Experience a revolutionary way to handle fat and optimize your clinical workflow with efficiency, safety and comfort with the new **Liposat® Pro plus**, specially designed for lipotransfer.

Key benefits at one glance

FAT Tubing specifically for fat transfer

Touch operation For TLA and fat infiltration

Connections for two foot switches

Minimum shear forces

Automated volume management

Noise optimisation

Delivery rate from 50 - 1000 ml/min

For small and large volume

Visualization of the infiltrated volume

Slims the workflow

Easy cleanability

Autostop technology



Surgest
Medical

Distribuidor Exclusivo en España
Teléfono de atención al cliente: 935 895 350
surgest.com

Five reasons

- 1** The liposuction handle offers the smoothest and most gentle way to remove fat.
- 2** Less effort is required to move the cannula through the tissue, ensuring a lower risk with maximum precision results.
- 3** Its anatomically shaped hand grip carries just the right balance point with a low weight, allowing the doctor to hold it comfortably for any length of time.
- 4** Work with the an innovative product with reliable and long-lasting performance.
- 5** The design allows for quick and easy cannula changing - providing an optimal, more efficient workflow.



▶ Our **Vibrasat® Pro** is the leading technology for vibration-assisted liposuctions. Its easy-to-use design eliminates fatigue for the doctor even during the longest procedures - while its technology maximizes patient results.

Key benefits at one glance

Boost function
Powerful and overheat protected motor Smooth running

Ergonomically cushioned **Anatomically shaped handle**
Reach up to 5,000 strokes/min with a stroke length of approx. 3mm

Convenient operation via foot switch possible QuickLock connection

Time saving



Surgest
Medical

Distribuidor Exclusivo en España
Teléfono de atención al cliente: 935 895 350
surgest.com

Four reasons

- 1** Keeps the temperature of the pre-heated infiltration fluid constant.
- 2** Integrates seamlessly into existing Möller Medical Aesthetics Solutions or operable independently. Tailored to your needs.
- 3** Keeps infiltration fluid warm during the procedure to maximize patient safety.
- 4** Its automatic temperature control simplifies workflows.



Hypothermia is a real risk for patients undergoing procedures where large amounts of fluid will be infiltrated into their bodies without a prior warming of the fluid. To reduce this risk Möller Medical offers you **Thermia Pro** - the **anti-hypothermia solution**. Thermia Pro ensures that the temperature of the fluid is maintained at a suitable temperature.

Key benefits at one glance

Built-in temperature
sensors

Simple to operate

Connections for two heating trays

Auto controlled

For safe TLA infiltration

Keeps the medium at body temperature

Stand-alone unit or mountable on the Workstations



Surgest
Medical

Distribuidor Exclusivo en España
Teléfono de atención al cliente: 935 895 350
surgest.com